

Great flexibility of use

The 9000 3D will soon cover the 10% of the cases that would require a complete view of the jaw at once in dentistry.


One can choose from the acquisition interface to conduct a detailed exam (76 μ m) locally or a global exam over a jaw (200 μ m). It all depends on the clinical indication.

We can make up to 3 exams of a given jaw and get them reconstructed into one, automatically...



Automatic acquisition protocol

Choice of 2 or 3 sites to be stitched together



Parameters:

- 200 μm reconstruction
- 23s x 3 = 69s scan time max
- 11s x 3 = 33s x-ray beam max
- reconstruction time < 2mn30s
- Outcome is a DICOM file

